

Guidelines for decontamination and declaration of contamination status

1. Ensure that the unit is disconnected from the mains supply.
2. Spray ALL external surfaces and tank with an appropriate disinfectant solution carefully following the manufacturers instructions.
3. Wipe all of the surfaces of the unit with a clean non-abrasive cloth and disinfectant solution and allow to dry.

Please complete the section below and overleaf and place inside the packaging of the unit you are returning.

Customer Name

Return Address

.....

.....

.....

Telephone

Email address

Customer ref

Model

Serial No

Fault

.....

Detail any damage to

the unit

...continued overleaf

Has the unit come into contact with any of the following contamination?

Blood, bodily fluids, biohazards, pathological
Items/samples, drugs, radioactive material,
hazardous chemicals?

YES / NO

Has the unit been decontaminated?

YES / NO

If YES please state with what product and how

.....
.....
.....

If NO please state why

.....
.....
.....

PLEASE NOTE. Walker Electronics Ltd reserves the right to charge a decontamination fee for units returned without this decontamination certificate or for a unit that has not been correctly decontaminated.

NOTE :

The unit should be returned whenever possible in its original packaging and within Walker Electronics Ltd accept no responsibility for damage in transit due to inadequate packaging.

The unit should be returned with the lid, basket and power cable. Any small items such as basket feet or lid pads will be replaced free of charge.

I confirm that the above is true and accurate to the best of my knowledge and that the unit has been prepared to ensure safe handling and transportation.

Name

Position

Signature

Date/...../.....