

RD5

General Purpose Cleaning Fluid



A concentrated biodegradable non-ammoniated general purpose cleaning fluid designed to clean, brighten and degrease most dental, podiatry and veterinary instruments.

Can also be used to remove polishing compounds, wax and general contamination from Jewellery and Watch Cases.

	1 Litre	5 Litre
Dilution rate	1 part concentrate to 9 parts of clean water	1 part concentrate to 9 parts of clean water
Makes a total of	10 litres	50 litres
Total weight including Container	1.08 kg	5.38 kg
Container dimensions W x D x H (mm)	125 x 70 x 165	190 x 140 x 285
MSDS link	https://walkerelectronics.co.uk/download/18/	